


AI221 ENGYNIOUS LIQUID LEVEL MONITOR APPLICATION DRAWING



FOR THIS APPLICATION
SET ITEM 60 TO 3

NOTE:
The Engynious level sensor provides
the current for input to the AI221

THIS IS AN UNCONTROLLED COPY - UNLESS STAMPED WITH CONTROLLED COPY STAMP

		JTL SYSTEMS LTD 41 KINGFISHER COURT, HAMBRIDGE RD, NEWBURY, BERKSHIRE: RG14 5SJ TEL: 01635 263646 FAX: 01635 263647	
ISSUE	DESCRIPTION	APP'D	DATE
1	ORIGINAL ISSUE	MJK	06/12
TITLE: <h2 style="text-align: center;">AI221</h2> <h3 style="text-align: center;">Liquid Level Application Drawing</h3>			
PROJ.	DRG No: 04003	SHT 1	OF 1
			ISS. 1